



Spencer Ligature Scissors - Straight/Blunt-Blunt/9cm



Used for cutting stitches or suture.

General Info

Product Number	14076-09
Category	Scissors & Micro Scissors
EAN	4061706001346



Find 14076-09 on finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	12 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	24 mm		
Cuttingedge	16 mm		
Length	9 cm		
Serrated	No		
Tips	Blunt/Blunt		
Tipshape	Straight		