



Lexer-Baby Scissors - Curved/Blunt-Blunt/10cm



General Info

Product Number	14079-10
Category	Scissors & Micro Scissors
EAN	4251909504052



Find 14079-10 on
finescience.com

Technical Info

Shipping Weight	20 g	Net Weight	17 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	24 mm		
Cuttingedge	16 mm		
Length	10 cm		
Serrated	No		
Tips	Blunt/Blunt		
Tipshape	Curved		