



Artery Scissors - Straight/Sharp-Ball Tip/11.5 cm



The ball tip prevents piercing of underlying tissue.

General Info

Product Number	14080-11
Category	Scissors & Micro Scissors > Fine Scissors
EAN	4251909504069



Find 14080-11 on
finescience.com

Technical Info

Shipping Weight	20 g	Net Weight	19 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	2		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	24 mm		
Cuttingedge	16 mm		
Length	11.5 cm		
Feature	1.5 mm Ball Diameter		
Serrated	No		
Tips	Sharp/Ball tip		
Tipshape	Straight		