



Dissecting Scissors - Heavy Blades/Angled to Side/Sharp-Sharp/10cm



Heavy blades for cutting fine bone.

General Info

Product Number	14082-09
Category	Scissors & Micro Scissors > Fine Scissors
EAN	4251909504083



Find 14082-09 on
finescience.com

Technical Info

Shipping Weight	20 g	Net Weight	17 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	12 mm		
Cuttingedge	8 mm		
Length	10 cm		
Serrated	No		
Tips	Sharp/Sharp		
Tipshape	Angled to side		