



A fine iris scissors, best suited for microdissection under magnification.

General Info

Product Number	14083-08
Category	Scissors & Micro Scissors > Fine Scissors
EAN	4251909504090



Find 14083-08 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	15 g	Net Weight	11 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Bladelength	15 mm		
Cuttingedge	10 mm		
Length	9 cm		
Serrated	No		
Tips	Sharp/Blunt		
Tipshape	Straight		