



Bonn Scissors - Curved/Sharp-Sharp/9cm



A fine iris scissors, best suited for microdissection under magnification.

General Info

Product Number	14085-08
Category	Scissors & Micro Scissors > Fine Scissors MedTech Line Most Popular > Scissors & Micro Scissors
EAN	4251909504120



Find 14085-08 on
finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	12 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Bladelength	15 mm		
Cuttingedge	10 mm		
Length	9 cm		
Serrated	No		
Tips	Sharp/Sharp		
Tipshape	Curved		