



Bonn Vascular Scissors - Straight/Sharp-Ball Tip/9cm



The ball tip prevents piercing of underlying tissue.

General Info

Product Number	14086-09
Category	Scissors & Micro Scissors > Fine Scissors
EAN	4251909504144



Find 14086-09 on
finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	13 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	15 mm		
Cuttingedge	10 mm		
Length	9 cm		
Feature	1.5 mm Ball Diameter		
Serrated	No		
Tips	Sharp/Ball tip		
Tipshape	Straight		