

F · S · T

Iris Scissors - Extra Narrow/Straight/Sharp-Sharp/10.5cm



Allgemein

Produktnummer	14088-10
Kategorie	Scissors & Micro Scissors > Fine Scissors
EAN	4251909504151



Finde 14088-10 auf
finescience.com

Technische Informationen

Versandgewicht	20 g	Nettogewicht	16 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Bladelength	24 mm		
Cuttingedge	16 mm		
Length	10.5 cm		
Serrated	No		
Tips	Sharp/Sharp		
Tipshape	Straight		