



Iris Scissors - Extra Narrow/Straight/Sharp-Sharp/10.5cm



General Info

Product Number	14088-10
Category	Scissors & Micro Scissors > Fine Scissors Most Popular > Scissors & Micro Scissors
EAN	4251909504151



Find 14088-10 on
finescience.com

Technical Info

Shipping Weight	20 g	Net Weight	16 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	24 mm		
Cuttingedge	16 mm		
Length	10.5 cm		
Serrated	No		
Tips	Sharp/Sharp		
Tipshape	Straight		