



Iris Scissors - Partly Polished/Straight/Sharp-Sharp/11cm



These scissors have undergone an extra vacuum hardening process that extends the life of their cutting edges, and they have an anti-reflection finish to reduce glare.

General Info

Product Number	14090-11
Category	Scissors & Micro Scissors > Fine Scissors Most Popular > Scissors & Micro Scissors
EAN	4251909504175



Find 14090-11 on
finescience.com

Technical Info

Shipping Weight	25 g	Net Weight	22 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Bladelength	24 mm		
Cuttingedge	16 mm		
Length	11 cm		
Serrated	No		
Tips	Sharp/Sharp		
Tipshape	Straight		