



Iris Scissors - Partly Polished/Curved/Sharp-Sharp/9cm



These scissors have undergone an extra vacuum hardening process that extends the life of their cutting edges, and they have an anti-reflection finish to reduce glare.

General Info

Product Number	14091-09
Category	Scissors & Micro Scissors > Fine Scissors Most Popular > Scissors & Micro Scissors
EAN	4251909504182



Find 14091-09 on
finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	12 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	24 mm		
Cuttingedge	16 mm		
Length	9 cm		
Serrated	No		
Tips	Sharp/Sharp		
Tipshape	Curved		