



Iris Scissors - Sandblasted/Curved/Sharp-Sharp/10.5cm



Made from martensitic stainless steel, a very hard form of steel crystalline, combined with molybdenum and vanadium. Molybdenum is an alloying agent used to form extremely hard steel with greater corrosion resistance, vanadium is the primary component of hardened carbon steel.

General Info

Product Number	14095-11
Category	Scissors & Micro Scissors > Fine Scissors
EAN	4251909504212



Find 14095-11 on
finescience.com

Technical Info

Shipping Weight	25 g	Net Weight	21 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Bladelength	24 mm		
Cuttingedge	16 mm		
Length	10.5 cm		
Serrated	No		
Tips	Sharp/Sharp		
Tipshape	Curved		