




Standard Surgical Scissors - Large Loops/Straight/Sharp-Blunt/14.5cm



These finger loops are approximately 30 larger than our standard surgical scissors.

General Info

Product Number	14101-14
Category	Scissors & Micro Scissors > Fine Scissors
EAN	4049356053991



Find 14101-14 on finescience.com

Technical Info

Shipping Weight	45 g	Net Weight	41 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	45 mm		
Cuttingedge	30 mm		
Length	14.5 cm		
Feature	Large loops		
Haslargeloops	Yes		
Serrated	No		
Tips	Sharp/Blunt		
Tipshape	Straight		