



Stevens Scissors - Adjustable Loops/Straight/Blunt-Blunt/9.5cm



The finger loops can be enlarged to fit almost any size.

General Info

Product Number	14108-09
Category	Scissors & Micro Scissors > Fine Scissors
EAN	4251909504243



Find 14108-09 on
finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	12 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	12 mm		
Cuttingedge	8 mm		
Length	10 cm		
Feature	Large loops		
Serrated	No		
Tips	Blunt/Blunt		
Tipshape	Straight		