



Mayo Scissors - ToughCut®/Straight/Blunt-Blunt/14.5cm



Allgemein

Produktnummer	14110-15
Kategorie	Scissors & Micro Scissors
EAN	4251909504250



Finde 14110-15 auf finescience.com

Technische Informationen

Versandgewicht	55 g	Nettogewicht	47 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	42 mm		
Cuttingedge	28 mm		
Length	14.5 cm		
Feature	ToughCut®		
Serrated	Yes		
Tips	Blunt/Blunt		
Tipshape	Straight		