



Mayo Scissors - ToughCut®/Straight/Blunt-Blunt/14.5cm



General Info

Product Number	14110-15
Category	Scissors & Micro Scissors
EAN	4251909504250



Find 14110-15 on finescience.com

Technical Info

Shipping Weight	55 g	Net Weight	47 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	42 mm		
Cuttingedge	28 mm		
Length	14.5 cm		
Feature	ToughCut®		
Serrated	Yes		
Tips	Blunt/Blunt		
Tipshape	Straight		