



Mayo Scissors - ToughCut®/Straight/Blunt-Blunt/17cm



Allgemein

Produktnummer	14110-17
Kategorie	Scissors & Micro Scissors
EAN	4251909504267



Finde 14110-17 auf finescience.com

Technische Informationen

Versandgewicht	70 g	Nettogewicht	67 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	48 mm		
Cuttingedge	32 mm		
Length	17 cm		
Feature	ToughCut®		
Serrated	Yes		
Tips	Blunt/Blunt		
Tipshape	Straight		