



Mayo Scissors - ToughCut®/Straight/Blunt-Blunt/17cm



General Info

Product Number	14110-17
Category	Scissors & Micro Scissors
EAN	4251909504267



Find 14110-17 on
finescience.com

Technical Info

Shipping Weight	70 g	Net Weight	67 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	48 mm		
Cuttingedge	32 mm		
Length	17 cm		
Feature	ToughCut®		
Serrated	Yes		
Tips	Blunt/Blunt		
Tipshape	Straight		