



Mayo Scissors - ToughCut®Curved/Blunt-Blunt/14.5cm



General Info

|                |                           |
|----------------|---------------------------|
| Product Number | 14111-15                  |
| Category       | Scissors & Micro Scissors |
| EAN            | 4251909504274             |



Find 14111-15 on [finescience.com](https://finescience.com)

Technical Info

|                        |                                 |                 |      |
|------------------------|---------------------------------|-----------------|------|
| Shipping Weight        | 50 g                            | Net Weight      | 48 g |
| Package Size           | 1 Pc.                           | Unit of Measure | EA   |
| Minimum Order Quantity | 1                               |                 |      |
| Alloymaterial          | Surgical Stainless Steel (Inox) |                 |      |
| Bladelength            | 42 mm                           |                 |      |
| Cuttingedge            | 28 mm                           |                 |      |
| Length                 | 14.5 cm                         |                 |      |
| Feature                | ToughCut®                       |                 |      |
| Serrated               | Yes                             |                 |      |
| Tips                   | Blunt/Blunt                     |                 |      |
| Tipshape               | Curved                          |                 |      |