



Metzenbaum Scissors - ToughCut®/Straight/Blunt-Blunt/14.5cm



General Info

Product Number	14116-14
Category	Scissors & Micro Scissors
EAN	4251909504298



Find 14116-14 on finescience.com

Technical Info

Shipping Weight	35 g	Net Weight	33 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	33 mm		
Cuttingedge	22 mm		
Length	14.5 cm		
Feature	ToughCut®		
Serrated	Yes		
Tips	Blunt/Blunt		
Tipshape	Straight		