



Metzenbaum Scissors - ToughCut®/Curved/Blunt-Blunt/14.5cm



Allgemein

Produktnummer	14117-14
Kategorie	Scissors & Micro Scissors
EAN	4251909504311



Finde 14117-14 auf finescience.com

Technische Informationen

Versandgewicht	35 g	Nettogewicht	35 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	33 mm		
Cuttingedge	22 mm		
Length	14.5 cm		
Feature	ToughCut®		
Serrated	Yes		
Tips	Blunt/Blunt		
Tipshape	Curved		