



Metzenbaum Scissors - ToughCut®/Curved/Blunt-Blunt/18cm



Allgemein

Produktnummer	14117-18
Kategorie	Scissors & Micro Scissors
EAN	4251909504328



Finde 14117-18 auf finescience.com

Technische Informationen

Versandgewicht	45 g	Nettogewicht	43 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	36 mm		
Cuttingedge	24 mm		
Length	18 cm		
Feature	ToughCut®		
Serrated	Yes		
Tips	Blunt/Blunt		
Tipshape	Curved		