



Metzenbaum Scissors - ToughCut®/Curved/Blunt-Blunt/18cm



General Info

Product Number	14117-18
Category	Scissors & Micro Scissors
EAN	4251909504328



Find 14117-18 on
finescience.com

Technical Info

Shipping Weight	45 g	Net Weight	43 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	36 mm		
Cuttingedge	24 mm		
Length	18 cm		
Feature	ToughCut®		
Serrated	Yes		
Tips	Blunt/Blunt		
Tipshape	Curved		