



Lister Scissors - ToughCut®/Angled to Side/Blunt-Blunt/14cm



Excellent skin cutting scissors. The longer, blunt blade prevents damaging tissue underneath skin without losing the cutting line.

General Info

Product Number	14131-14
Category	Scissors & Micro Scissors
EAN	4251909504342



Find 14131-14 on finescience.com

Technical Info

Shipping Weight	50 g	Net Weight	47 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	39 mm		
Cuttingedge	26 mm		
Length	14 cm		
Feature	ToughCut®		
Serrated	Yes		
Tips	Blunt/Blunt		
Tipshape	Angled to side		