



Iris Scissors - Straight / Sharp-Sharp/11.5cm/Left-Handed



Allgemein

Produktnummer	14160-10
Kategorie	Scissors & Micro Scissors > Fine Scissors
EAN	4251909504359



Finde 14160-10 auf finescience.com

Technische Informationen

Versandgewicht	20 g	Nettogewicht	18 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	27 mm		
Cuttingedge	18 mm		
Length	11.5 cm		
Islefthanded	Yes		
Serrated	No		
Tips	Sharp/Sharp		
Tipshape	Straight		