



Iris Scissors - Straight / Sharp-Sharp/11.5cm/Left-Handed



General Info

Product Number	14160-10
Category	Scissors & Micro Scissors > Fine Scissors
EAN	4251909504359



Find 14160-10 on
[finescience.com](https://www.finescience.com)

Technical Info

Shipping Weight	20 g	Net Weight	18 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	27 mm		
Cuttingedge	18 mm		
Length	11.5 cm		
Islefthanded	Yes		
Serrated	No		
Tips	Sharp/Sharp		
Tipshape	Straight		