



Iris Scissors - Curved/Sharp-Sharp/11.5cm/Left-Handed



General Info

Product Number	14161-10
Category	Scissors & Micro Scissors > Fine Scissors
EAN	4251909504366



Find 14161-10 on
finescience.com

Technical Info

Shipping Weight	20 g	Net Weight	17 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	27 mm		
Cuttingedge	18 mm		
Length	11.5 cm		
Islefthanded	Yes		
Serrated	No		
Tips	Sharp/Sharp		
Tipshape	Curved		