



Kaye Scissors - Adjustable Loops/Curved/Blunt-Blunt/10.5cm



These scissors have semi-sharp outer edges to better facilitate dissection.

General Info

Product Number	14171-10
Category	Scissors & Micro Scissors
EAN	4251909504373



Find 14171-10 on finescience.com

Technical Info

Shipping Weight	20 g	Net Weight	19 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	21 mm		
Cuttingedge	14 mm		
Length	10.5 cm		
Haslargeloops	Yes		
Serrated	Yes		
Tips	Blunt/Blunt		
Tipshape	Curved		