



One edge is serrated to minimize slippage of tissue during cutting. Curved handle helps when cutting from an awkward angle.

General Info

Product Number	14173-12
Category	Scissors & Micro Scissors
EAN	4251909504380



Find 14173-12 on
finescience.com

Technical Info

Shipping Weight	20 g	Net Weight	17 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	18 mm		
Cuttingedge	12 mm		
Length	11.5 cm		
Serrated	Yes		
Tips	Sharp/Sharp		
Tipshape	Curved		