



Bonn Scissors - Highly Polished/Straight/Sharp-Sharp/9cm



Allgemein

Produktnummer	14184-09
Kategorie	Scissors & Micro Scissors > Fine Scissors
EAN	4251909504397



Finde 14184-09 auf finescience.com

Technische Informationen

Versandgewicht	15 g	Nettogewicht	11 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	15 mm		
Cuttingedge	10 mm		
Length	9 cm		
Tips	Sharp/Sharp		
Tipshape	Straight		