



Standard Surgical Scissors - Tungsten Carbide/Straight/Sharp-Blunt/14.5cm



General Info

Product Number	14501-14
Category	Scissors & Micro Scissors > Surgical Scissors
EAN	4251909504441



Find 14501-14 on finescience.com

Technical Info

Shipping Weight	45 g	Net Weight	42 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox) with tungsten carbide		
Bladelength	51 mm		
Bladematerial	Tungsten Carbide		
Cuttingedge	34 mm		
Length	14.5 cm		
Feature	Tungsten Carbide		
Serrated	No		
Tips	Sharp/Blunt		
Tipshape	Straight		