



Mayo Scissors - Tungsten Carbide/Straight/Blunt-Blunt/17cm



General Info

Product Number	14510-17
Category	Scissors & Micro Scissors
EAN	4251909504519



Find 14510-17 on
finescience.com

Technical Info

Shipping Weight	70 g	Net Weight	67 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox) with tungsten carbide		
Bladelength	48 mm		
Bladematerial	Tungsten Carbide		
Cuttingedge	32 mm		
Length	17 cm		
Feature	Tungsten Carbide		
Serrated	No		
Tips	Blunt/Blunt		
Tipshape	Straight		