



Mayo Scissors - Tungsten Carbide/Curved/Blunt-Blunt/14.5cm



General Info

Product Number	14511-15
Category	Scissors & Micro Scissors
EAN	4251909504526



Find 14511-15 on finescience.com

Technical Info

Shipping Weight	50 g	Net Weight	47 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox) with tungsten carbide		
Bladelength	42 mm		
Bladematerial	Tungsten Carbide		
Cuttingedge	28 mm		
Length	15 cm		
Feature	Tungsten Carbide		
Serrated	No		
Tips	Blunt/Blunt		
Tipshape	Curved		