



Mayo Scissors - Tungsten Carbide & ToughCut®/Straight/Blunt-Blunt/14.5cm



General Info

Product Number	14512-15
Category	Scissors & Micro Scissors
EAN	4251909504540



Find 14512-15 on finescience.com

Technical Info

Shipping Weight	55 g	Net Weight	48 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox) with tungsten carbide		
Bladelength	42 mm		
Bladematerial	Tungsten Carbide		
Cuttingedge	28 mm		
Length	15 cm		
Feature	Tungsten Carbide & ToughCut®		
Serrated	Yes		
Tips	Blunt/Blunt		
Tipshape	Straight		