




Metzenbaum Fino Scissors - Tungsten Carbide/Curved/Blunt-Blunt/14.5cm



General Info

Product Number	14519-14
Category	Scissors & Micro Scissors
EAN	4251909504588



Technical Info

Shipping Weight	25 g	Net Weight	25 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox) with tungsten carbide		
Autoclavesafe	Yes		
Bladelength	33 mm		
Bladematerial	Tungsten Carbide		
Cuttingedge	22 mm		
Length	14.5 cm		
Feature	Tungsten Carbide		
Serrated	No		
Tips	Blunt/Blunt		
Tipshape	Curved		