



Metzenbaum Fino Scissors - Tungsten Carbide/Curved/Blunt-Blunt/18cm



General Info

Product Number	14519-18
Category	Scissors & Micro Scissors
EAN	4251909504595



Find 14519-18 on finescience.com

Technical Info

Shipping Weight	35 g	Net Weight	35 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox) with tungsten carbide		
Autoclavesafe	Yes		
Bladelength	36 mm		
Bladematerial	Tungsten Carbide		
Cuttingedge	24 mm		
Length	18 cm		
Feature	Tungsten Carbide		
Serrated	No		
Tips	Blunt/Blunt		
Tipshape	Curved		