



Iris Scissors - Tungsten Carbide/ToughCut®/Straight/Sharp-Sharp/11.5cm



MOST
POPULAR



TUNGSTEN CARBIDE
TOUGH CUT



SERRATED

General Info

Product Number	14558-11
Category	Scissors & Micro Scissors > Fine Scissors MedTech Line Most Popular > Scissors & Micro Scissors
EAN	4251909504625



Find 14558-11 on
finescience.com

Technical Info

Shipping Weight	25 g	Net Weight	24 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox) with tungsten carbide		
Bladelength	24 mm		
Bladematerial	Tungsten Carbide		
Cuttingedge	16 mm		
Length	11.5 cm		
Feature	Tungsten Carbide & ToughCut®		
Serrated	Yes		
Tips	Sharp/Sharp		
Tipshape	Straight		