

F · S · T

Iris Scissors - Tungsten Carbide/ToughCut®/Curved/Sharp-Sharp/9cm



General Info

Product Number	14559-09
Category	Scissors & Micro Scissors > Fine Scissors MedTech Line
EAN	4251909504632



Find 14559-09 on
finescience.com

Technical Info

Shipping Weight	20 g	Net Weight	13 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox) with tungsten carbide		
Bladlength	24 mm		
Bladematerial	Tungsten Carbide		
Cuttingedge	16 mm		
Length	9 cm		
Feature	Tungsten Carbide & ToughCut®		
Serrated	Yes		
Tips	Sharp/Sharp		
Tipshape	Curved		