



Iris Scissors - Tungsten Carbide/Straight/Sharp-Sharp/11.5cm



General Info

Product Number	14568-12
Category	Scissors & Micro Scissors > Fine Scissors MedTech Line Most Popular > Scissors & Micro Scissors
EAN	4251909504663



Find 14568-12 on finescience.com

Technical Info

Shipping Weight	25 g	Net Weight	24 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox) with tungsten carbide		
Bladelength	24 mm		
Bladematerial	Tungsten Carbide		
Cuttingedge	16 mm		
Length	11.5 cm		
Feature	Tungsten Carbide		
Serrated	No		
Tips	Sharp/Sharp		
Tipshape	Straight		