




Metzenbaum Baby Scissors - Tungsten Carbide/Straight/Blunt-Blunt/11.5cm



General Info

Product Number	14574-11
Category	Scissors & Micro Scissors MedTech Line
EAN	4251909504700



Find 14574-11 on finescience.com

Technical Info

Shipping Weight	25 g	Net Weight	23 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox) with tungsten carbide		
Bladelength	30 mm		
Bladematerial	Tungsten Carbide		
Cuttingedge	20 mm		
Length	11 cm		
Feature	Tungsten Carbide		
Serrated	No		
Tips	Blunt/Blunt		
Tipshape	Straight		