



Metzenbaum Mini Scissors - Tungsten Carbide/Curved/Blunt-Blunt/9cm



General Info

Product Number	14575-09
Category	Scissors & Micro Scissors MedTech Line
EAN	4251909504717



Find 14575-09 on
finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	12 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox) with tungsten carbide		
Bladelength	18 mm		
Bladematerial	Tungsten Carbide		
Cuttingedge	12 mm		
Length	9 cm		
Feature	Tungsten Carbide		
Serrated	No		
Tips	Blunt/Blunt		
Tipshape	Curved		