



Metzenbaum Baby Scissors - Tungsten Carbide/Curved/Blunt-Blunt/11.5cm



General Info

|                |                                           |
|----------------|-------------------------------------------|
| Product Number | 14575-11                                  |
| Category       | Scissors & Micro Scissors<br>MedTech Line |
| EAN            | 4251909504724                             |



Find 14575-11 on [finescience.com](https://finescience.com)

Technical Info

|                        |                                                       |                 |      |
|------------------------|-------------------------------------------------------|-----------------|------|
| Shipping Weight        | 25 g                                                  | Net Weight      | 22 g |
| Package Size           | 1 Pc.                                                 | Unit of Measure | EA   |
| Minimum Order Quantity | 1                                                     |                 |      |
| Alloymaterial          | Surgical Stainless Steel (Inox) with tungsten carbide |                 |      |
| Bladelength            | 30 mm                                                 |                 |      |
| Bladematerial          | Tungsten Carbide                                      |                 |      |
| Cuttingedge            | 20 mm                                                 |                 |      |
| Length                 | 11 cm                                                 |                 |      |
| Feature                | Tungsten Carbide                                      |                 |      |
| Serrated               | No                                                    |                 |      |
| Tips                   | Blunt/Blunt                                           |                 |      |
| Tipshape               | Curved                                                |                 |      |