



Bonn-Gawenda Scissors - ToughCut®/Straight/Sharp-Sharp/10.5cm



Introducing our innovative Bonn-Gawenda Scissors from the Exclusive Line, blending precision cutting with elongated handles to produce seamless microscopy work. Our ToughCut® system features micro serrated cutting edges ensuring an exact cut without slippage, and new larger finger rings which enhance control and comfort compared to other standard iris scissors.

General Info

Product Number	14854-10
Category	Exclusive Line Scissors & Micro Scissors > Fine Scissors
EAN	4251909511982



Find 14854-10 on finescience.com

Technical Info

Shipping Weight	73 g	Net Weight	19.5 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Bladelength	15 mm		
Cuttingedge	10 mm		
Length	10.5 cm		
Feature	ToughCut®		
Haslargeloops	Yes		
Tips	Sharp/Sharp		
Tipshape	Straight		