



Iris Scissors - CeramaCut®/Straight/Sharp-Sharp/9cm



Our CeramaCut® scissors combine ToughCut® technology with Tungsten Carbide edges and a ceramic coating for one of the best, cleanest cutting scissors available on the market today. One of the hardest materials on Earth, ceramic is used in surgical instruments to form edges that last longer than stainless steel and even longer than Tungsten Carbide. Ceramic also has a high slide capacity which reduces surface tension and results in a much cleaner cut than metal instruments.

General Info

Product Number	14958-09
Category	Scissors & Micro Scissors > Fine Scissors MedTech Line Most Popular > Scissors & Micro Scissors
EAN	4251909504731



Find 14958-09 on finescience.com

Technical Info

Shipping Weight	15 g	Net Weight	11 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Ceramic coated Surgical Stainless Steel (Inox) with tungsten carbide		
Autoclavesafe	Yes		
Bladelength	24 mm		
Bladematerial	Tungsten Carbide		
Cuttingedge	16 mm		
Length	9 cm		
Feature	CeramaCut®		
Serrated	Yes		
Tips	Sharp/Sharp		
Tipshape	Straight		