



Cohan-Vannas Spring Scissors - Curved/Sharp/6cm/6mm Cutting Edge



Allgemein

Produktnummer	15000-01
Kategorie	Scissors & Micro Scissors > Micro Scissors
EAN	4251909504786



Finde 15000-01 auf finescience.com

Technische Informationen

Versandgewicht	25 g	Nettogewicht	4 g
Paketgröße	1 Pc.	Maßeinheit	EA
Mindestbestellmenge	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cuttingedge	6 mm		
Length	6 cm		
Serrated	No		
Tipdiameter	0.1 mm		
Tips	Sharp/Sharp		
Tipshape	Curved		