



Cohan-Vannas Spring Scissors - Curved/Sharp/6cm/6mm Cutting Edge



Long cutting blades for a more controlled cut.

General Info

Product Number	15000-01
Category	Scissors & Micro Scissors > Micro Scissors
EAN	4251909504786



Find 15000-01 on finescience.com

Technical Info

Shipping Weight	25 g	Net Weight	4 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cuttingedge	6 mm		
Length	6 cm		
Serrated	No		
Tipdiameter	0.1 mm		
Tips	Sharp/Sharp		
Tipshape	Curved		