



Vannas Spring Scissors - Straight/Sharp/8cm/2mm Cutting Edge



These scissors have been designed to maximize visibility under the microscope by eliminating unnecessary bulk in the blades. These Vannas spring scissors have the smallest blades ever produced.

General Info

Product Number	15000-03
Category	Scissors & Micro Scissors > Micro Scissors
EAN	4251909504809



Find 15000-03 on finescience.com

Technical Info

Shipping Weight	25 g	Net Weight	4 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cuttingedge	2 mm		
Length	8 cm		
Serrated	No		
Tipdiameter	0.05 mm		
Tips	Sharp/Sharp		
Tipshape	Straight		