




Vannas Spring Scissors - Straight/Sharp/8cm/2.5mm Cutting Edge



These scissors have been designed to maximize visibility under the microscope by eliminating unnecessary bulk in the blades. A short taper makes these scissors ideal for small, precise cuts.

General Info

Product Number	15000-08
Category	Scissors & Micro Scissors > Micro Scissors Most Popular > Scissors & Micro Scissors
EAN	4251909504823



Find 15000-08 on
finescience.com

Technical Info

Shipping Weight	25 g	Net Weight	4 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cuttingedge	2.5 mm		
Length	8 cm		
Serrated	No		
Tipdiameter	0.05 mm		
Tips	Sharp/Sharp		
Tipshape	Straight		