



Cohan-Vannas Spring Scissors - Curved/Blunt/6cm/6mm Cutting Edge



Long cutting blades for a more controlled cut.

General Info

Product Number	15000-11
Category	Scissors & Micro Scissors > Micro Scissors
EAN	4251909504847



Find 15000-11 on
finescience.com

Technical Info

Shipping Weight	25 g	Net Weight	4 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Autoclavesafe	Yes		
Cuttingedge	6 mm		
Length	6 cm		
Serrated	No		
Tipdiameter	0.3 mm		
Tips	Blunt/Blunt		
Tipshape	Curved		