



Vannas-Tübingen Spring Scissors - Straight/Sharp/8.5cm/5mm Cutting Edge



General Info

Product Number	15003-08
Category	Scissors & Micro Scissors > Micro Scissors Most Popular > Scissors & Micro Scissors
EAN	4251909504885



Find 15003-08 on finescience.com

Technical Info

Shipping Weight	25 g	Net Weight	4 g
Package Size	1 Pc.	Unit of Measure	EA
Minimum Order Quantity	1		
Alloymaterial	Surgical Stainless Steel (Inox)		
Cuttingedge	5 mm		
Length	8.5 cm		
Serrated	No		
Tipdiameter	0.05 mm		
Tips	Sharp/Sharp		
Tipshape	Straight		